

"Keep Honey Fraud at Bay"

Honey Authenticity Testing by $^1\text{H-NMR}$



Léa Heintz, Product Manager
Bruker BioSpin GmbH
09.09.2019, Montréal, Canada



Honey-Profiling overview



Honey-Profiling is an application module of the **FoodScreener platform**



- Standardized procedures for sample preparation and measurement to ensure data reproducibility and the usage of common data analysis tools.
- Sample measurement and data analysis are fully automatized
- Sample preparation requires 2,5 mg of honey diluted by 1 to 4 in water followed by pH adjustment

Sugar syrups in honey



- According to EU-Directive of 2001 and Codex Alimentarius, honey should be exempt of foreign sugars.
- Foreign sugars can be:
 - a deliberate addition to honey in order to stretch the product.
 - a consequence of **bee feeding** with sugar syrups.
- Feeding bees with industrial sugar syrups: \approx 40 times less expensive than feeding with honey.
- Sugar syrups need to be handled with care and be given outside of production periods.
- Famille Michaud, the first reseller of honey in France tests all the honeys they buy with their NMR instrument. **They find sugar syrups in 10% of the French honeys tested.**


https://www.francetvinfo.fr/replay-jt/france-2/20-heures/video-les-abeilles-nourries-au-sirop-de-betterave_3521387.html



Detection of Sugar Syrups with specific methods



- Foreign enzymes
 - Particular enzyme
- Single Markers
 - Rice syrup
 - Beet syrup



With the courtesy of QSI GmbH

- AOAC EA-IRMS
 - C4-type sugars

- Detection of **25%** of adulterated honeys in Honey-Profiling Database
- ¹³C EA-IRMS only: acceptance of approximately 3/4 adulterated honeys.

➤ You need to know what you are looking for

➤ Tests are easy to deceive



Detection of Sugar Syrups with NMR Honey-Profiling



- **Quantitative analysis** of many markers: up to 60 markers per sample
- The values are compared to the reference thresholds
- Reference thresholds have been determined and validated thanks to a large database of authentic and known adulterated honeys from various countries and botanical varieties:
 - 16100 authentic samples / 1900 known adulterated samples
 - > 50 countries & 100 botanical varieties.
 - Monofloral & polyfloral honeys from a single country source.
 - Blends of honeys from different countries.
 - Industrial honeys (or “baker honey”)

Why using a database ?

Validation of markers and related thresholds



Example of Fructose/Glucose ratio

- Independently of the analytical technique used (NMR, MS, HPLC-UV...) and of the parameter measured (unknown compound or identified molecule like e.g. mannose):


To define what is "normal" it is mandatory to have a set of representative samples, to make sure markers chosen are not naturally present in a specific type of honey.

➤ AVOID THE RISK OF FALSE POSITIVE

- That is even more important when looking to low concentrated compounds (like e.g. secondary metabolites of plants)
- In NMR Honey-Profiling, the markers are all in the spectral regions related to sugars, which is less sensitive to origin.

Detection of false declaration of origin

Pollen analysis



- Identification and counting pollen grains by microscopy
- Requires a trained expert and is time consuming
- Not working on pollen-filtered honeys.
- Not able to differentiate the real pollens from exogenous pollen grains.

Detection of false declaration of origin

Marker compounds



- Compounds commonly used as markers can often be synthetically produced and added into the honey.
- E.g. DHA and MGO are cheap chemicals
- Non-compliance of marker compounds can be used to detect a fraud
- However, **compliance of 1 or few markers are not sufficient to prove authenticity**

Compound	Value	Unit	LOQ	Reference Range	Flag
3-phenyllactic acid	682	mg/kg	300	388  2300	●
dihydroxyacetone (DHA)	378	mg/kg	20	235  3513	●
methylglyoxal (MGO)	243	mg/kg	30	110  1486	●
					-

Riddle of how 1,700 tons of manuka honey are made... but 10,000 are sold

23 Aug, 2016 1:30pm

3 minutes to read




Research by the main honey producers' organisation in New Zealand has revealed that 1,700 tonnes of manuka are produced there each year. Photo / Brett Phibbs

Detection of false declaration of origin NMR Honey-Profiling



- Based on the **complete chemical composition of the honey** observed by 1H-NMR
- Statistical analysis of the NMR spectrum: untargeted buckets/variables.
- Differentiation of the supposed/declared origin or variety with all other ones present in the Database (50 countries / 100 varieties).
- Statistical models have been validated by Monte Carlo cross validation (Criteria: TP > 98%)
 - ***Very difficult/ impossible to deceive***
 - ***No need for an expert and fast***
 - ***Applicable to pollen-filtered honey***




Verification of Labeling of Origin on pollen-filtered honey




^1H -Spectra of USA/ARG. (40:60) Honey: raw vs filtered

p8711_10_1
p8746_10_1



p8711_10_1
p8746_10_1



Verification of origin possible even on pollen-filtered honey with NMR !

Additional information with NMR

Honey-Profiling



- **Detection of atypical samples**
 - Statistical comparison of the NMR profile with “normal” profiles from the same floral source
 - Potential to detect new frauds at early stage.
- **Honey constituent's quantification**
 - Sugars, organic acids, amino acids, freshness and quality criteria
 - Absolute quantification
 - Comparison to reference values (obtained by our Database samples)
- **Regulated parameters in EU-directive / Codex Alimentarius**
 - HMF, glucose+fructose, sucrose
 - Absolute quantification
 - Conclusion according to directive.

Why using NMR and Honey-Profiling



- An **all-in-one method**, comprehensive authenticity test:
 - Detection of sugar syrups
 - Verification of labelling/declaration of origin (country and botanical source)
 - Detection of atypical samples
 - Quality control (composition, freshness, regulated parameters)
- **deeply validated method**, relying on a huge Database of 18000 well-characterized samples
- Several labs (including Bruker BAS lab) accredited **ISO17025** for Honey-Profiling
- no NMR expertise required. Easy to operate, fully automated, no regular instrument cleaning
- Fast: 20 – 25 min / sample
- Difficult to deceive
- Wide and ww network of partners and users

NMR Honey-Profilng network



- A **worldwide network** of partners and customers.
- Honey-Profilng is a **centralized database** to which states official laboratories can contribute. The goal is to have a database of authentic material covering all its natural variations and to define **common/harmonized purity criteria**
- Possibility to develop new models according to country-specific needs **much quicker**, and with **better performance**.

i.e to characterize one country it is necessary to have samples from other countries to know what are the differences that should be looked for.

- **Round tables** have been kicked-off to discuss the NMR purity criteria defined amongst honey experts

